Lab Exercise #6

rev 1

Objectives

- Learn about the Serial Peripheral Interface (SPI)
- Implement a demonstration SPI loop-back program to observe SPI operation
- Write routines to initialize and then gather data from an STM LIS2MDL magnetometer

References

- [1] NXP Kinetis KL25 processor sub-family data sheet (pdf) on class web page.
- [2] NXP Kinetis KL25 processor sub-family reference manual (pdf) on class web page.
- LIS2MDL Magnetometer data sheet (pdf) on class web page

Lab 6 Problem Statement

Implement the SPI Loopback test routines listed on pages 224 to 226 in the textbook.

Design Flow

Lab Part A

- 1) Create a new folder and project for SPI loopback testing.
- 2) Copy files from prior project(s), as needed, to the new folder
- 3) Implement the loopback test as shown in the textbook Monitor the SPI signals with three scope channels.
- 4) Submit your loopback project to D2L as Lab6a

Lab Part B

- 1) Create a new folder and project for SPI communication with the LIS2MDL magnetometer
- 2) Create a program to read data from two axes of the magnetometer and display it in a terminal window at an update rate of about twice a second.
- 3) Submit lab part B files to D2L as Lab6b

To Turn In

- In the "comment header" of your main.c files report success, failure, or other observations
- Submit your lab6a main.c file to a D2L drop box for lab6a.
- Zip up your complete lab6a uVision project and submit to the lab6a D2L drop box
- Submit your lab6b main.c file to a D2L drop box for lab6b.
- Zip up your complete lab6b uVision project and submit to the lab6b D2L drop box